

Public Health England Microbiology Services Porton Down Salisbury Wiltshire SP4 OJG

## **Certificate of Containment Testing**

Containment Testing of Rotor S-4x750 (5895 120.105-00) with Plate Buckets (5895 124.119-00\*) and Caps (5820 743.301-00\*) in an Eppendorf Bench Top Centrifuge

Report No. 14/043 A

Report Prepared For: Eppendorf AG, Hamburg, Germany

**Issue Date:** 17<sup>th</sup> February 2015

## **Test Summary**

Rotor S-4x750 (5895 120.105-00) with Plate Buckets (5895 124.119-00\*) and Caps (5820 743.301-00\*) was containment tested in an Eppendorf bench top centrifuge, using Annex AA of IEC 61010-2-020:2006 (2<sup>nd</sup> Ed.). The sealed buckets were shown to contain a spill.

Report Written By

Name: Ms Anna Moy

Title: Biosafety Scientist

Anna Way

Report Authorised By

Name: Mrs Sara Speight

Title: Senior Biosafety Scientist

Please be aware that the use of the Royal Coat of Arms is highly restricted and cannot be copied. Please do not put the PHE logo on your website or use our name to endorse your products. Any reference to PHE needs to be approved by us before it can be used.

<sup>\*</sup> Part no. will form part of catalogue number 5895 128.009; 5895 124.003; 5895 125.000

<sup>&</sup>lt;sup>#</sup> Part no. will form part of catalogue number 5820 748.001