

**Compatibility of Combitips advanced with standard laboratory tubes**

Combitips advanced	Eppendorf Safe-Lock Tubes				Conical tubes			Eppendorf Deepwell Plates			
	0.5 mL	1.5 mL	2.0 mL	5.0 mL	15 mL	25 mL	50 mL	96 wells/ 500 µL	96 wells/ 1000 µL	96 wells/ 2000 µL	384 wells/ 200 µL
0.1 mL	+	+	+	+	-	-	-	+	+	+	+
0.2 mL	+	+	+	+	-	-	-	+	+	+	+
0.5 mL	+	+	+	+	-	-	-	+	+	+	-
1 mL	+	+	+	+	-	[+]	-	+	+	+	-
2.5 mL	[+]	[+]	[+]	[+]	[+]	[+]	[+]	+	[+]	[+]	-
5 mL	+	+	+	+	[+]	[+]	[+]	+	+	+	-
10 mL	+	+	+	-	-	+	[+]	+	+	+	-
25 mL	-	+	+	-	-	-	-	+	+	+	-
50 mL	-	+	+	-	-	-	-	+	+	+	-

[+] = improved compatibility compared to the Combitip plus generation